



E-xpand[®]

Stentgraft Balloon Catheter

- Designed to assist in stentgraft procedures e.g. E-vita thoracic or E-vita abdominal
- Compliant balloon
- Radiopaque markers for exact placement
- Short inflation and deflation times
- Compatible with 0.035" (0,89 mm) guide wire
- Minimum sheath diameter 12F (4 mm)

Stentgraft Balloon Catheter

Ordering information

Quick and reliable delivery of standard products in Europe within 24 hours.

Call our central office: +49 (0) 7471 / 922-0

E-xpand Stentgraft Balloon Catheter

ø Balloon (mm)	ø Shaft (F/mm)	Usable Length (cm)	Minimum Sheath Diameter (F/mm)	Catalogue No.
10-50	8/2.7	100	12/4	85XX0050N35-00

