



E[®]-asy
plus

Introducer Sheath

- Eases the introduction of endovascular catheter and stentgraft systems into the vascular system
- Simplest single-hand control
- Continuously adjustable hemostatic valve guarantees full range sealing from 0F to 24F
- Coiled tubing for maximum kink resistance and flexibility
- Hydrophilic coating enables safe and secure insertion
- Available in sizes of 16F, 18F, 20F, 22F and 24F



Introducer Sheath

Ordering information

Quick and reliable delivery in Europe within 24 hours.

Call our central office: +49(0)7471/922-0

E-asy plus Introducer Sheath

Length (cm)	Sheath size (F/mm)	Catalogue No.
20	16/5.3	74XX2016N-35
20	18/6.0	74XX2018N-35
20	20/6.7	74XX2020N-35
20	22/7.3	74XX2022N-35
30	20/6.7	74XX3020N-35
30	22/7.3	74XX3022N-35
30	24/8.0	74XX3024N-35



JOTEC AG
Luzernerstr. 91
CH-5630 Muri
T +41 (0) 56 670 91 - 00
F +41 (0) 56 670 91 - 02
www.jotec.net

JOTEC GmbH
Lotzenäcker 23
D-72379 Hechingen
T +49 (0) 74 71 / 922 - 0
F +49 (0) 74 71 / 922 - 100

JOTEC s.r.l.
Vicolo Valtellina n.15
I-20092 Cinisello Balsamo (Mi)
T +39 02 30310061
F +39 02 66014455

JOTEC CARDIOVASCULAR SL
Av. Castilla 2, P. E. San Fernando
Edificio Italia – planta baja
E-28830 San Fernando de Henares, Madrid
T +34 (0) 91 - 648 70 20
F +34 (0) 91 - 656 17 73

