



E[®]-wire

Guide Wire

- Coated, extra stiff 0.035" guide wire
- Short and highly flexible J-tip with outstanding X-ray visibility
- Available in two different lengths: 190 and 300 cm
- For use in any catheter based diagnostic and interventional procedures

Guide Wire

Ordering information

Quick and reliable delivery in Europe within 24 hours.

Call our central office: +49 (0) 7471 / 922-0

E-wire Guide Wire

Length (cm)	Length flexible tip (cm)	Diameter ("/mm)	Catalogue No.
190	6	0.035/0.89	76XX1935N-06
300	6	0.035/0.89	76XX3035N-06

3 units/box