



E-njoy[®]

Self Expandable Peripheral Stent

- Laser-cut Nitinol stent
- Diameters from 5 to 10 mm, lengths from 30 to 100 mm
- Flexible stent design with minimal foreshortening allowing effective vessel reconstruction
- 6 Tantal markers for good X-ray visibility and exact placement
- 6F delivery system for use with 0.035" (0.89 mm) guide wires
- Two different working lengths: 80 and 130 cm
- Easy pull-back delivery system allowing frictionless and controlled stent positioning
- Special catheter material (braided tube) ensures high flexibility and precise torquability
- Atraumatic tip for smooth vessel crossing

Self Expandable Peripheral Stent

Ordering information

Quick and reliable delivery in Europe within 24 hours.

Call our central office: +49 (0) 7471 / 922-0

E-njoy Self Expandable Peripheral Stent

Catheter Length (cm)	Diameter (mm)	Stent Length (mm)	Catalogue No.
80	5	30	81PS3005
80	5	40	81PS4005
80	5	60	81PS6005
80	6	40	81PS4006
80	6	60	81PS6006
80	6	80	81PS8006
80	7	40	81PS4007
80	7	60	81PS6007
80	7	80	81PS8007
80	8	40	81PS4008
80	8	60	81PS6008
80	8	80	81PS8008
80	9	40	81PS4009
80	9	60	81PS6009
80	9	80	81PS8009
80	10	40	81PS4010
80	10	60	81PS6010
80	10	80	81PS8010
80	10	100	81PS0010
130	5	30	81PL3005
130	5	40	81PL4005
130	5	60	81PL6005
130	6	40	81PL4006
130	6	60	81PL6006
130	6	80	81PL8006
130	7	40	81PL4007
130	7	60	81PL6007
130	7	80	81PL8007
130	8	40	81PL4008
130	8	60	81PL6008
130	8	80	81PL8008
130	9	40	81PL4009
130	9	60	81PL6009
130	9	80	81PL8009
130	10	40	81PL4010
130	10	60	81PL6010
130	10	80	81PL8010
130	10	100	81PL0010

Product Specifications

E-njoy Self Expandable Stent

Recommended introducer sheath	6F (2 mm)
Maximum guide wire	0.035" (0.89 mm)